

LA-UR-15-28284

Approved for public release; distribution is unlimited.

Title: Introduction to neutron spectrometers

Author(s): Mclean, Thomas Donaldson

Intended for: Neutron training class at CSU (Fort Collins, CO) Oct27-29 2015

Issued: 2015-10-22

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the Los Alamos National Security, LLC for the National Nuclear Security Administration of the U.S. Department of Energy under contract DE-AC52-06NA25396. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.



Introduction to neutron spectrometers

Tom McLean, LANL

CSU neutron detection class
Fort Collins, CO
Oct. 27-29 2014

UNCLASSIFIED

Introduction: talk outline

- Discussion (brief) of neutron spectrometers for neutrons < 20 MeV
 - Bonner spheres
 - ROSPEC
 - Scintillator-based instruments

UNCLASSIFIED

Slide 2

Introduction


- Spectroscopic methods of neutron detection provide:
 - Neutron spectral data
 - Neutron dose by folding in fluence-to-dose conversion coefficients
- However spectrometers are more complicated than rem meters
 - Therefore more expensive and usually less reliable
 - More complicated to operate
 - Typically require sophisticated real-time pulse shape discrimination
 - Usually, a method of energy calibration is required
 - Typically they are not portable
 - Typically do not provide real-time information as off-line unfolding or deconvolution of the pulse height data is required.
 - requiring a well-characterized instrument response function

UNCLASSIFIED

Slide 3

Bonner spheres


- Developed in 1960's
- Based on measurements taken with PE spheres of different diameters (3"-12")
- Typically uses a ^3He counter but also TLDs and foils for use in pulsed fields.
- Time consuming and labor intensive measurement process
- Data is unfolded to give the neutron spectrum
 - Spectrum consists of 53 energy bins based on only 10 or fewer measurements !!!



UNCLASSIFIED

Slide 4

Neutron Spectrometers : BSS



Data also applies to the response of a PE-moderated rem meter

UNCLASSIFIED


Neutron Spectrometers : BSS

- Broad response functions
 - Calculated using Monte Carlo techniques (e.g. MCNP)
- Possible scenarios
 - Neutron escapes after one or more interactions
 - Neutron is thermalized but captured on H ($H(n,\gamma)D$)
 - Neutron is thermalized and enters counter
- Increasing PE sphere diameter :
 - Allows thermalization of increasingly higher energy neutrons
 - Decrease in detection efficiency for lower energy neutrons
 - 8"-10" diameter response approx. dose equivalent curve

UNCLASSIFIED

Neutron Spectrometers : BSS

Typical BSS pulse height spectrum



Uses pulse height to discriminate against gammas

Integrate count rate above "A" to get response of the sphere in the neutron field being measured.


UNCLASSIFIED

Neutron Spectrometers : BSS




- Unfolding the BSS data is not straightforward due to the under determined nature of the measurements (i.e. fewer measurements than neutron energy bins)
- Therefore a unique solution is not possible
- However, a solution spectrum that is most consistent with an initial guess spectrum is possible
- A variety of unfolding codes have been developed over the years
- Notable new codes include MAXED and GRAVEL developed by PTB

UNCLASSIFIED

Neutron Spectrometers : BSS




Count rate data



UNCLASSIFIED


Neutron Spectrometers : BSS



MAXED
solution

UNCLASSIFIED

Neutron Spectrometers : BSS




MAXED
solution

UNCLASSIFIED

Neutron Spectrometers : BSS

Bonner Sphere spectra

SRS data : Bare Cf spectrum dependence on distance



UNCLASSIFIED

Neutron Spectrometers : ROSPEC



6 gas proportional counters on rotating platform

- 4 based on hydrogenous detector medium (proton recoil)
 - 3 2 inch diameter (different pressures)
 - 1 6 inch diameter
- 2 ^3He counters
 - 2 2 inch diameter
 - Bare for thermal neutrons
 - ^{10}B -shielded for epi-thermal detection

UNCLASSIFIED

Neutron Spectrometers : ROSPEC

- Pulse height spectrum
 - 256 channels per counter
 - Each counter has limited dynamic range
 - Gamma contribution at low energies
 - Wall effects at higher energies

- Unfold (deconvolute)
 - Using Benjamin code (SPEC) 1970 vintage
 - Adequate for our purposes (improvements have been made)
 - Better codes available
 - Recoil proton response functions
 - Ideally rectangular but resolution and wall effects perturb

UNCLASSIFIED


Neutron Spectrometers : ROSPEC

CTR#	CTR name	Diameter (inches)	Gas	Pressure (atm)	Range (MeV)
0	SP2-1 ^a	2	H ₂	0.75	0.05-0.25
1	SP2-4	2	H ₂	4.0	0.15-0.70
2	SP2-10	2	H ₂	10.0	0.40-1.50
3	SP6	6	P10	5.0	1.20-5.0
4		2	³ He	0.08	<0.5x10 ⁻⁶
5		2	³ He	0.5	0.5x10 ⁻⁶ -1x10 ⁻²

notes: a - Spherical 2 inch - 1 atmosphere b - 90% Ar + 10% CH₄

UNCLASSIFIED

Neutron Spectrometers : ROSPEC




Monoenergetic neutron

$$E_n = 1.2 \text{ MeV}$$

UNCLASSIFIED


Neutron Spectrometers : ROSPEC

ROSPEC II NEUTRON SPECTRUM



UNCLASSIFIED

Neutron Spectrometers : ROSPEC




Monoenergetic
neutron

$$E_n = 146 \text{ keV}$$

UNCLASSIFIED


Neutron Spectrometers : ROSPEC



D₂O-Cf

UNCLASSIFIED

Neutron Spectrometers : ROSPEC



UNCLASSIFIED

Neutron Spectrometers : ROSPEC


- Not necessary to get “publication” quality spectrum to determine dose accurately
- Because :
 - Even unsophisticated unfolding codes conserve fluence
 - Dose curves also insensitive to fluence distribution above 1MeV

UNCLASSIFIED

Neutron Spectrometers : ROSPEC

Realistic field spectra

Realistic field @ 150 cm $E_p=3\text{MeV}$ May 6 2004




UNCLASSIFIED

Neutron Spectrometers : BSS vs ROSPEC

Comparison of ROSPEC and BSS spectra

LANL NFIA : D₂O-moderated Cf at 200cm




UNCLASSIFIED

Neutron Spectrometers : BSS vs ROSPEC

Comparison of ROSPEC and BSS spectra

LANL NFIA : D₂O-moderated Cf at 200cm



UNCLASSIFIED

Neutron Spectrometers : N-Probe



- Manufactured by BTI
- Detectors :
 - 1"x1" NE213 (E>800 keV)
 - ^3He (boron covered)
- Portable (incl. analyzer unit)
- Active (n, γ) discrimination (>1000:1 at low rates)
- Real time dose data
- Crude (but adequate) spectroscopy


UNCLASSIFIED

Neutron Spectrometers : N-Probe

- NE213 (aka BC501A)
 - Xylene-based liquid scintillator
 - Proton recoil and reactions on Carbon
 - Two principal decay times
 - 30ns & 300ns
 - Slow component fraction is proportional to particle dE/dx
 - Excellent PSD properties
 - Non-linear light output
 - dE/dX dependent saturation effect
 - Limited lifetime due to O_2 quenching

UNCLASSIFIED

Neutron Spectrometers : N-Probe




UNCLASSIFIED

Neutron Spectrometers : N-Probe

Pulse shape discrimination


$E_n=19.6$ MeV



UNCLASSIFIED


Neutron Spectrometers : N-Probe

N-PROBE : LOW GAIN NE213 SPECTRUM



UNCLASSIFIED


Neutron Spectrometers : N-Probe



UNCLASSIFIED

Neutron Spectrometers : N-Probe

Comparison of N-Probe and ROSPEC



UNCLASSIFIED

Neutron Spectrometers : SSS

- Manufactured by BTI

- Detector :

 - 16 3x3mm plastic scintillators


 - (recoil protons produced)

- Portable (incl. analyzer)

- High-energy extension to ROSPEC

- Passive n, γ discrimination

- Low efficiency




Glass light pipe

PMT

UNCLASSIFIED


Neutron Spectrometers : SSS

SSS MICROSPEC-2 DATA




UNCLASSIFIED

Neutron Spectrometers : SSS



UNCLASSIFIED

Neutron Spectrometers : SSS



UNCLASSIFIED

Summary of spectrometer properties

Instrument	Detectors	Energy range (MeV)	Max. dose rate (mrem/hr) for bare Cf	Comments
ROSPEC	Gas proportional + $^3\text{He}(n,p)$	Thermal – 5 MeV	200	Pulse height n,γ discrimination
Bonner spheres	$^3\text{He}(n,p)$	Thermal – 30 MeV	1000	Pulse height n,γ discrimination
N-Probe	Liquid scintillator + $^3\text{He}(n,p)$	Thermal – 20 MeV	100	Pulse shape n,γ discrimination. Real-time dose
SSS	Plastic scintillator	5 - 15 MeV	>1000	Pulse height n,γ discrimination

UNCLASSIFIED

Other methods of interest

- $\text{Cs}_2\text{LiYCl}_6:\text{Ce}$ (CLYC) scintillator
 - Gamma spectroscopy (high resolution)
 - Thermal neutron detection
 - Fast neutron detection
- LLNL-developed scintillators commercialized by Eljen Technology
 - Excellent n/gamma discrimination
 - Better than NE213
- Time of flight techniques
 - Take advantage of relatively slow neutron velocities

UNCLASSIFIED

Slide 37